

## WD6010

---



10 DIN BASKETS DOUBLE DOOR WASHER DISINFECTOR  
DOUBLE DOOR WASHER DISINFECTOR FOR CSSD  
[Download PDF](#)

---

- ✓ washing up to 95 °C - antibacterial thermal disinfection
- ✓ capacity of 12 DIN baskets
- ✓ programs: 20 default + 20 custom
- ✓ multi-language TFT colour touch screen on both sides
- ✓ functioning protected by different levels password
- ✓ integrated built-in panel printer (optional)
- ✓ USB port for downloading washing cycles log
- ✓ Ethernet port for LAN connection (optional)
- ✓ n. 2 drop down full glass doors
- ✓ doors interlocked for complete separation dirty/clean side
- ✓ n. 2 spray arms, one on the top and on the bottom of the chamber
- ✓ high pressure hot air drying system with absolute 99.999% HEPA filters
- ✓ high efficiency Ecoslim steam condenser (optional)
- ✓ washing rack recognition system (optional)
- ✓ peristaltic pump for liquid detergent and level sensor
- ✓ peristaltic pump for liquid neutralizer and level sensor
- ✓ possibility of installing up to 2 additional peristaltic pumps
- ✓ detergent cabinet with removable tray for storing up to 3 x 5 liters cans + 1 x 1 liter can
- ✓ washing chamber made of AISI 316L stainless steel
- ✓ external panels made of AISI 304 stainless steel
- ✓ external dimensions LxDxH: 650 x 720 x 1950 mm
- ✓ internal working dimensions LxDxH: 550 x 620 x 685 mm
- ✓ cold water, hot water and demineralized water connections
- ✓ floor drain connection
- ✓ drain pump (optional)

- ✓ electrical connection: 3/N/PE 400 V - 50 Hz - 13 kW max
- ✓ 60 Hz version available
- ✓ without racks

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350  
Phone +39 0522 82111 - fax +39 0522 821453 - [service.instruments@smeg.it](mailto:service.instruments@smeg.it)

