

WD5090TFO



WASHER DISINFECTOR WITH DRYING SYSTEM, EN15883 APPROVED - FULL OPTIONAL
PROFESSIONAL LINE - WASHER DISINFECTOR 90 CM WITH DRYING SYSTEM
[Download PDF](#)

- ✓ washing up to 95 °C - antibacterial thermal disinfection
- ✓ programs: 20 default + 10 custom
- ✓ n. 3 microprocessors programmable electronic control
- ✓ multi-language graphic display, 128 x 64 pixel
- ✓ serial port RS232 for connection to external printer or PC
- ✓ integrated printer
- ✓ forced hot air drying system with absolute 99.999% HEPA filter
- ✓ high efficiency steam condenser Ecoslim, with cold water injection
- ✓ peristaltic pump for liquid detergent and level sensor
- ✓ peristaltic pump for liquid neutralizer and level sensor
- ✓ peristaltic pump for liquid disinfectant and level sensor
- ✓ peristaltic pump for liquid lubricant and level sensor
- ✓ redundant dosage system control with flow meters
- ✓ washing chamber and internal door in AISI 316L stainless steel
- ✓ external panels made of AISI 304 stainless steel
- ✓ detergent cabinet made of AISI 304 stainless steel, 4 x 5 liters tanks capacity
- ✓ upper level can holder basket
- ✓ external dimensions LxDxH: 900 x 630 x 850 mm
- ✓ internal working dimensions LxDxH: 525 x 490 x 570 mm
- ✓ incorporated water softener with salt sensor
- ✓ cold water, hot water and demineralized water connections
- ✓ electronic door lock system for user safety
- ✓ electrical connection: 3/N/PE 400 V - 50 Hz - 7 kW max (switchable to single-phase 230 V)
- ✓ 60 Hz version available
- ✓ without trolleys

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350
Phone +39 0522 8211 - fax +39 0522 821453 - service.instruments@smeg.it

