

Haier Biomedical

Intelligent Protection of Life Science

Combined Refrigerator and Freezer

HYCD-282A



Refrigerator & Freezer Combination

Independent display for refrigerator section and freezer section



Alarm Functions

High/low temperature, power failure, sensor error, low battery, high/low ambient temperature, and remote alarm



Standard Stainless Steel Interior

Both refrigerator and freezer chamber are stainless steel interior, durable, easy to clean and anti-corrosion

Ergonomic Design



Lockable design



Intelligent alarm



Foot/Caster design



Test hole design

Specifications

HYCD-282A

Technical Data

| | |
|---------------------|---|
| Type | Combined Refrigerator and Freezer |
| Cabinet Type | Upright |
| Climate Class | 4 |
| Cooling Type | Refrigerator: Forced air cooling Freezer: Direct cooling |
| Defrost Mode | Refrigerator: Auto Freezer: Manual |
| Refrigerant | CFC-Free |
| Noise Level (dB(A)) | 43 |

Performance

| | |
|-----------------------|---------------------------------------|
| Temperature Range(°C) | Refrigerator: 2~8 Freezer: -20~-40 |
|-----------------------|---------------------------------------|

Control

| | |
|------------|----------------|
| Controller | Microprocessor |
| Display | LED |

Electrical Data

| | |
|-----------------------|------------|
| Power Supply(V/Hz) | 220~240/50 |
| Power(W) | 400 |
| Electrical Current(A) | 2.8 |

Dimensions

| | | |
|-------------------------------|---|--|
| Capacity(L/Cu.Ft) | Refrigerator:185/6.5 | Freezer:97/3.4 |
| Net/Gross Weight(approx) | 145/160(kg) 319.7/352.7(lbs) | |
| Interior Dimensions(W*D*H) | Refrigerator:605*510*720(mm) Freezer:515*465*440(mm) | 23.8*20.0*28.3(in) 20.3*18.3*17.3(in) |
| Exterior Dimensions(W*D*H) | 736*660*1810(mm) | 29.0*26.0*71.3(in) |
| Packing Dimensions(W*D*H) | 775*720*1930(mm) | 30.5*28.3*76(in) |
| Container load (20'/40'/40'H) | 21/45/45 | |



Functions

| | |
|----------------------|---|
| High/Low Temperature | Y |
| Remote Alarm | Y |
| Power Failure | Y |
| Sensor Error | Y |
| Low Battery | Y |
| High Ambient Temp | Y |

Accessories

| | |
|-----------------|-----------------------------------|
| Caster | Y |
| Foot | Y |
| Test Hole | Y |
| Shelves/Drawers | Refrigerator: 3/1 Freezer: 0/2 |
| USB Interface | Y |

Others

| | |
|-------------|----|
| Certificate | CE |
|-------------|----|

Haier Biomedical UK

Address: 2 Birch House, Brambleside, Bellbrook
Industrial Estate, Uckfield, East Sussex,
TN22 1QQ, United Kingdom
TEL: +44(0) 1825 744 696 (Head Office)
www.haierbiomedical.co.uk



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Scan 2D barcode with your mobile device to learn more information about Haier Biomedical